

Basic Water Still




Unique Design & Safety Feature
with hydraulic temperature controller and a water float switch provides automatic operation.

Protect heater from over heating by cut-out when water is low level

Easy to Clean Up
by removing the condenser from the body

All Stainless Steel with Double Walled Chamber
finished with powder coated two colors



WD Series Basic Water Still

- ❖ Space Saving Compact and Practical Design
- ❖ Maximum 7.5 l/hour Production Capacity
- ❖ Economical Energy and Water Consumption



Model	WD-1004	WD-1008	
	WD-2004F	WD-2008F	
Capacity	3.5 l / hour	7.5 l / hour	
Dimensions (mm)	380 (W) x 290 (D) x 580 (H)	380 (W) x 365 (D) x 640 (H)	
Heater	3 kW	6 kW	
Cooling Water	Approx. 30 ~ 50 l	Approx. 60 ~ 80 l	
Material	Inner	Stainless Steel	
	Outer	Stainless Steel /w Powder Coating	
	Heater	Stainless Steel Pipe	
Controller	Hydraulic Thermostat		
Electrical Requirements	220V, 50/60 Hz		
Switch	Main Power with Pilot Lamp and Water Float Switch		
Sealing	Silicone Gasket		
Safety	Low Water & Over Heating Cut-Out		

Order Number	Descriptions	
WD-1004	Basic Water Still	3.5 l /hour
WD-2004F†	Basic Water Still /w Prefilter	3.5 l /hour
WD-1008	Basic Water Still	7.5 l /hour
WD-2008F†	Basic Water Still /w Prefilter	7.5 l /hour



WD-2004F Basic Water Still /w Prefilter

☒ Custom made order is available
Specifications can be changed without prior notice for quality improvement.