



WD-1008 & WD-1004 Basic Water Still

Unique Design & Safety Feature

with hydraulic temperature controller and a water float switch provides automatic operation.

Protect heater from over heating by cut-out when water is low level

Easy to Clean Up

by removing the condenser from the body

All Stainless Steel with Double Walled Chamber

finished with powder coated two colors

WD Series Basic Water Still

- ❖ Space Saving Compact and Practical Design
- ❖ Maximum 7.5 l/hour Production Capacity
- ❖ Economical Energy and Water Consumption

Our New Technology will provide you with the Water Still with Air cooled condensing.

No waste of Water

Coming Soon!!



Model	WD-1004	WD-1008
	WD-2004F	WD-2008F
Capacity	3.5 l / hour	7.5 l / hour
Dimensions (mm)	380 (W) x 290 (D) x 580 (H)	380 (W) x 365 (D) x 640 (H)
Heater	3 kW	6 kW
Cooling Water	Approx. 30 ~ 50 l	Approx. 60 ~ 80 l
Material	Inner Outer Heater	Stainless Steel Stainless Steel /w Powder Coating Stainless Steel Pipe
Controller	Hydraulic Thermostat	
Electrical Requirements	220V, 50/60 Hz	
Switch	Main Power with Pilot Lamp and Water Float Switch	
Sealing	Silicone Gasket	
Safety	Low Water & Over Heating Cut-Out	

Order Number	Descriptions
WD-1004	Basic Water Still 3.5 l /hour
WD-2004F+	Basic Water Still /w Prefilter 3.5 l /hour
WD-1008	Basic Water Still 7.5 l /hour
WD-2008F+	Basic Water Still /w Prefilter 7.5 l /hour



WD-2004F Basic Water Still /w Prefilter

☒ Custom made order is available
Specifications can be changed without prior notice for quality improvement.