


# Vertical Laminar Flow Clean Bench



## **Specification :**

Overall Dimension : 1300 x 750 x2000 mm

Working Dimension : 1200 x 645 x 610 mm

Display : Digital Air Flow

Pre-Filter : Recyclable,3-30  $\mu\text{m}$  Partical Removal

Main Filter : HEPA (High Efficiency Particulate Air) Filter 0.3  $\mu\text{m}$  Particle Removal 99.997%

Air Velocity : 0.3-0.5 m/s

Illumination : FL Lamp 28 W (750-800 Lumen)

Sterilization : UV Lamp 30 W

Material Worktable : Stainless Steel 304

Outer : Powder Coated Steel

Door : Tempered Safety Glass

Electric Supply : 220 Volt 50 Hz

Certification : CE, ISO 9001, ISO 14001, ISO 13485

